



Prepare plate with 95 samples

- Create plate in Geneious/Biocode LIMS
- Enter sample IDs by scanning 2D barcodes
- Print plate layout with sample positions



scan



Make notes on plate layout sheet:

- 1st individuum will become morphological voucher
- 2nd individuum from same location for Biobank

Also note any other special requirements.

Print labels for morphological vouchers (in DiversityCollection)

Accession: ZFMK-TS-250426 Germany, Schorndorfen, Euchleben, 51.881287°N/10.891703°E, Thüringen, Euchleben, im Gemischtbauch 18.VI.2013 leg. Kopetz, Andreas	<i>Ceryphilia marginella</i> Motsch., 1858 det. Kopetz, Andreas	Accession: ZFMK-TS-250427 Germany, Schorndorfen, Euchleben, 51.881287°N/10.891703°E, Thüringen, Euchleben, Wipfhaus NW, Altm 17.VI.2013 leg. Kopetz, Andreas	<i>Bombus terrestris</i> (Goebl., 1793) det. Kopetz, Andreas
Accession: ZFMK-TS-250430 Germany, 2304 NN, Erfurt, 51.0242°N/10.0226°E, Thüringen, Erfurt, Roter Berg, Zooпарк 30.VI.2013 leg. Kopetz, Andreas	<i>Berytus inceptor</i> (Goebl., 1863) det. Kopetz, Andreas	Accession: ZFMK-TS-250433 Germany, 2304 NN, Erfurt, 51.0241°N/10.0226°E, Thüringen, Erfurt, Roter Berg, Zooпарк 30.VI.2013 leg. Kopetz, Andreas	<i>Stenobothrus cylindricollis</i> (Fahrs., 1840) det. Hartmann, Matthias
Accession: ZFMK-TS-250441 Germany, 2804 NN, Nordhausen, 51.511215°N/10.478213°E, Thüringen, Salzgerath, Schellenberg, Südfang 4.VI.2013 leg. Kopetz, Andreas	<i>Acrobium nitidum</i> F., 1792 det. Kopetz, Andreas	Accession: ZFMK-TS-250442 Germany, 2804 NN, Nordhausen, 51.511215°N/10.478213°E, Thüringen, Salzgerath, Schellenberg, Südfang 4.VI.2013 leg. Kopetz, Andreas	<i>Ceryphilia marginella</i> Motsch., 1858 det. Kopetz, Andreas

Photograph samples



Fill S-block with lysis buffer and ProteinaseK



Transfer tissue of
each specimen into
corresponding well

(make notes on
printout)

Lab space for subsampling



Lab space for subsampling

Microscope & light source



Lab space for subsampling



Fresh gloves

Lab space for subsampling



S-block filled with ATL buffer and Proteinase K and silicon pad


Lab space for subsampling



2 forceps (1 soft and 1 sharp)
and petri dish

tissues & ethanol for cleaning

Lab space for subsampling




Lab space for subsampling

lighter & burner



Lab space for subsampling



Box with specimens, corresponding printout of plate, pen and additional box for already processed samples


Lab space for subsampling



Additional fresh vials for the voucher
and empty box for Biobank subsamples

Lab space for subsampling

paper labels (DC) for vials with voucher and scissors



Lab space for subsampling



Garbage bin

Larger arthropods: remove up to 3 legs from same (right) side of body




Larger arthropods: remove up to 3 legs from same (right) side of body

1 leg in S-block for lysis
(chopped up for a better inflow of Proteinase K)



Larger arthropods: remove up to 3 legs from same (right) side of body


1 leg in S-block for lysis
(chopped up for a better inflow of Proteinase K)



1 leg in labelled vial for Biobank

Larger arthropods: remove up to 3 legs from same (right) side of body

1 leg in S-block for lysis
(chopped up for a better inflow of Proteinase K)



1 leg in labelled vial for Biobank



Specimen with paper label in new vial as voucher (collection)

Examples of different sizes



→ 1 leg for lysis, 1 leg for Biobank



→ 2 legs for lysis, 1 leg for Biobank



→ 3 legs for lysis, no leg for Biobank



→ Whole body lysis

Do not take the pedipalps instead of legs!



→ 1 leg for lysis, 1 leg for Biobank



→ 2 legs for lysis, 2 legs for Biobank




→ 4 legs for lysis, no leg for Biobank



→ Whole body lysis


Tissue samples from larger animals



Piece of tissue (e.g. muscle)
sent from collector

Other tissue sources: toe or quill (reptiles and birds)
tip of tail (reptiles and rodents)
fin (fish)

Cover S-block with silicon pad and put into lysis overnight:
gently shaking at 300rpm and 56°C, ~14-16 hours




Next day: recover small individuals that have been in lysis completely → voucher vial with paper label

Fill 5 additional S-blocks and elution plate with buffers
(pipetting by hand or using the „EpMotion 5075“ robot)



Wash 1: 650µl Buffer AW1
Wash 2: 500µl Buffer AW1
Wash 3: 500µl Buffer AW2
Wash 4: 500µl Buffer AW2
Wash 5: 500µl H₂O
Elution: 200µl Buffer AE



In a 50ml tube prepare magnetic bead solution:

- 22ml Buffer AL
- 22ml Isopropanol
- 3.2 ml MagAttract (vortex very well before using)

→ Vortex solution, then add 450 μ l to each well

Load all blocks into BioSprint 96 machine and start according program.



Prepare PCR immediately after extraction, before storing the DNA

Using the Qiagen Multiplex Kit prepare a mix for 100 reactions:

- 440 μ l H₂O
- 200 μ l Q-Solution
- 1000 μ l Multiplex-Mix
- 80 μ l forward-Primer: LCO1490-JJ (5'-CHACWAAYCATAAAGATATYGG-3')
- 80 μ l reverse-Primer: HCO2198-JJ (5'-AWACTTCVGGRTGVCCAAARAATCA-3')

→ 1800 μ l Mix in total → fill 18 μ l in each well of a PCR-plate

Then transfer 2 μ l DNA of each well from the elution plate into the corresponding well of the PCR-plate (using a multi-channel pipette)

Cover plate with slide and put directly into PCR machine



Transfer remaining DNA from elution plate to long term storage system (FluidX)

Used PCR program is a combination of a 'touchdown' and a 'step-up' routine:

hot start Taq activation:	15 min at 95 °C
first cycle set (15 repeats):	35 s denaturation at 94°C
	90 s annealing at 55°C (-1°C per cycle)
	90 s extension at 72°C
second cycle set (25 repeats):	35 s denaturation at 94°C
	90 s annealing at 40°C
	90 s extension at 72°C
final elongation	10 min at 72 °C
cooling	12°C (∞)

Unpurified PCR products are sent for bidirectional Sanger sequencing to BGI (Hong Kong)

What we have in the end:



Vial with the specimen voucher



Temporary freezer until barcode has been validated



Permanent morphological collection

What we have in the end:



Vial with the specimen voucher



Temporary freezer until barcode has been validated



Permanent morphological collection

Subsample (leg or whole animal) for Biobank



What we have in the end:



Vial with the specimen voucher



Temporary freezer until barcode has been validated



Permanent morphological collection

Subsample (leg or whole animal) for Biobank



DNA (in storage tubes) → send to the different partners

